

mumu PLUS Extra


NITRILE EXAMINATION AND PROTECTIVE GLOVES



Latex-Free
Tactile Sensitivity
Chemical Resistance
Protection Against Viruses

Mumu Plus Extra examination gloves are made of high quality nitrile, which reduces the risk of allergies as they do not contain natural latex and powder.

Resistance against Bacteria and Fungi - **PASS**
Resistance against Virus - **PASS**



SPECIFICATION

Glove Details	Specification
Material	Nitrile (NBR)
Protein Level	Latex free&Powder free
Surface	Finger Textured
Shape	Ambidextrous
Colour	Blue
Sterilization	Non Sterile
Shelf Life	5 Years
Cuff Finished	Beaded

PHYSICAL PROPERTIES

Glove Details	Properties
Force at Break – Before & After Aging (Newtons)	min. 6
Elongation at Break – After Aging (%)	min. 400
Tensile Strength – Before & After Aging (MPA)	min. 14
Freedom from Holes	AQL 1.5
Length (mm)	Min. 240
Cuff Thickness (mm)	0.07 +/- 0,02
Palm Thickness (mm)	0.07 +/- 0,02
Finger Thickness (mm)	0.09 +/- 0,02

SIZING

Size	Palm Width
Extra Small	76 ± 3 mm
Small	84 ± 3 mm
Medium	94 ± 3 mm
Large	105 ± 3 mm
Extra Large	113 ± 3 mm

QUANTITIES

Inner Box Dimension	Carton Dimension	20 FT	40 HC
200x110x60 mm	455x315x210 mm	960 Cartons	2236 Cartons
Packaging: 100 pcs per box, 20 inner boxes carton			
Inner Box Dimension	Outer Carton	20 FT	40 HC
220x120x70 mm	365x250x230 mm	1400 Cartons	3220 Cartons
Packaging : 200 pcs per box, 10 inner boxes per carton			

Product Standards: The product is tested and manufactured to the following standards: EN455-1, EN455-2, EN455-3, EN455-4. EN ISO 374-1, EN 374-2, EN 16523-1+A1, EN 374-4, EN ISO 374-5, EN ISO 10993-1, EN 420, ASTM D6978-05.

Storage: Store in a cool dry place. Avoid excessive heat. Opened box should be shielded from exposure to direct sun or fluorescent lighting.

Disposal: Dispose of gloves as clinical waste. Transit case and box can be recycled as paper or disposed of as clinical waste. Please follow local regulation and guidance.

Shelf Life: Five(5) years from date of manufacture.

Country of origin: Malaysia

Notice: After use wearers should visually check the gloves and remove any contamination from the outer surface before removing the gloves from the hand. Alternatively carefully peel the gloves off the hand so that the contaminated glove outer does not touch your skin.

Registering authority: In Europe the gloves are CE marked (the notified body SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, DIS YN2P, Ireland, number 2777) indicating compliance with EU Regulation 2016/425.